



Product Discription

Mortuary tables are used for autopsies, embalming, dressing the body for a funeral, temporary body storage, and for teaching purposes. These tables are specifically designed to keep embalming and bodily fluids away from the mortician or medical professional. A well-constructed autopsy table is critical to maintaining a safe and clean work environment within your mortuary or funeral home.

GA2003FA Forensic Autopsy Table Features

1. Arc seamless AISI 304 stainless steel tabletop: The tabletop adopts an arc seamless connection design to avoid the accumulation of bacteria and dirt, easy to clean and disinfect, and ensure hygiene.
2. Laser scale screen isolation board: The tabletop is equipped with a laser scale screen isolation board, which is convenient for accurate measurement and recording, ensuring the professionalism and standardization of the anatomical operation.
3. Hair isolation system: A hair isolator is installed in the sink to effectively prevent hair and small debris from entering the drain pipe, avoid pipe blockage, and keep the drainage system smooth.
4. Hot and cold water adjustment and intelligent induction faucet: Equipped with hot and cold water adjustment function and induction faucet, it can be connected to the water heater, support automatic manual switching, ensure the water temperature is suitable, and easy to operate.
5. Mobile spray gun: Equipped with a mobile high-pressure spray gun, the water pressure is adjustable, which can quickly and effectively wash the dirt on the surface of the dissection table to meet high-standard cleaning needs.
6. Table lifting function: The tabletop has a 150mm (180mm optional) lifting stroke adjustment function, which is convenient for anatomical operations at different angles and helps to observe the details of inconvenient areas.
7. Multi-point cleaning head design: 8 cleaning heads are installed on the side of the table, which can quickly remove stains and keep the table clean more efficiently during the cleaning process.
8. Tissue crushing system: A tissue crusher is installed under the sink to crush the residual tissue after dissection and discharge it smoothly through the drainage pipe to ensure the cleanliness and sterility of the operation area.
9. Optional mobile tray: A mobile tray can be installed above the table to temporarily place the internal organs removed during the dissection process, which is convenient for operators to organize and observe (optional function).
10. Negative pressure ventilation system: Equipped with a negative pressure ventilation system, the upper part is ventilated through the vents on both sides of the table, and the openings around the base are used for lower ventilation, which can timely extract harmful gases generated during the dissection process and improve the air quality of the operating environment.
11. Multi-angle rotation function: The table can rotate 45° forward and backward, which is suitable for multi-angle shooting and teaching, and enhances the flexibility of teaching practice and observation.
12. Safe power protection: The side of the countertop is equipped with a leakage protection device and a 220V, 50Hz waterproof socket to ensure the safety of electricity during use.

Options: headrest, instrument trolley.

GA2003FA Forensic Autopsy Table Specifications:

Size: 2600* 800* 850mm

Table lifting stroke 150 (180)mm

Rotation 45° forward and backward

